

## Two Year Limited Warranty

HoMedics, Inc., guarantees this product free from defects in material and workmanship for a period of two years from the date of purchase, except as noted below.

This HoMedics product warranty does not cover damage caused by misuse or abuse; accidents; the attachment of any unauthorized accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of HoMedics. HoMedics shall not be responsible for any type of incidental, consequential or special damages. All implied warranties, including but not limited to those implied warranties of fitness and merchantability, are limited in the total duration of two years from the original purchase date.

To obtain warranty service on your Foot Pro Ultra product either hand deliver or mail the unit and your dated sales receipt (as proof of purchase), postpaid, along with a check or money order in the amount of \$5.00 payable to HoMedics, Inc., to cover handling.

Upon receipt, HoMedics will repair or replace, as appropriate, your Foot Pro Ultra product and return it to you, postpaid. Warranty is solely through HoMedics Service Center. Service of this product by anyone other than HoMedics Service Center voids warranty.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

**Mail To:**  
**HoMedics**  
**Service Center Dept. 168**  
**3000 Pontiac Trail**  
**Commerce Township MI**  
**48390**

**Email:**  
**cservice@homedics.com**

©2000 HoMedics, Inc. or its affiliated companies, all rights reserved.  
HoMedics® is a registered trademark of HoMedics Inc. or its affiliated companies.  
Foot Pro™ is a trademark of HoMedics Inc. or its affiliated companies.  
Z0255  
IB-AK3-A

**2**year  
limited warranty

## FootPro™ Ultra Luxury Foot Massager



## Instruction Manual and Warranty Information

AK-3

**2**year  
limited warranty

## IMPORTANT SAFETY INSTRUCTIONS:

When using an electrical appliance, basic precautions should always be followed, including the following:  
**READ ALL INSTRUCTIONS BEFORE USING**

### **DANGER** - To reduce the risk of electric shock:

- Always unplug the appliance from the electrical outlet immediately after using and before cleaning.
- Do not reach for an appliance that has fallen into water. Unplug it immediately.
- Do not use while bathing or showering.
- Do not place or store appliance where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.

### **WARNING** - To reduce the risk of burns, fire, electric shock or injury to persons:

- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts or attachments.
- Close supervision is necessary when this appliance is used by, on or near children, invalids or disabled persons.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by HoMedics; specifically any attachments not provided with the unit.
- Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return it to HoMedics Service Center for examination or repair.

**Caution:  
In case of  
pregnancy,  
diabetes or  
illness, consult  
your doctor  
before using  
massager.**

- Keep cord away from heated surfaces.
- Never block the air openings of the appliance or place it on a soft surface, such as a bed or couch, where the air openings may be blocked. Keep air openings free of lint, hair, etc.
- Never use while sleeping or fall asleep while using the massager.
- Never drop or insert any object into any opening.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electrical shock or injury to persons.
- Do not carry this appliance by supply cord or use cord as handle.
- To disconnect, depress all control buttons to the "OFF" position, then remove plug from outlet.
- This appliance is designed for household use only.
- Use heated surfaces carefully. May cause serious burns. Do not use on sensitive areas.

## SAVE THESE INSTRUCTIONS

**Caution – Please read all instructions carefully before operating.**

- If you have any concerns regarding your health, consult your doctor before using the Foot Pro Ultra.
- Individuals with pacemakers should not use this product.
- If you experience pain in a muscle or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.
- Never use Foot Pro Ultra directly on open wounds or discolored areas.

- Do not use for more than 15 minutes at a time. Excessive use could lead the product to overheat. Should this occur, discontinue use and allow the unit to cool before operating.
- If heat function is left on for a prolonged period of time, lightly touch bulb surface to ensure that heat level is not great. If too great, turn heat function off, allowing surface to cool.
- Never stand or bear heavy weight on the unit. Unit must be used while seated only. Standing or bearing heavy weight upon the unit will cause irreparable damage.
- Never try to open or stick objects into any part of the Foot Pro Ultra. Damages of this nature are not covered under the HoMedics warranty.
- The massage unit must be in the "OFF" position before plugging into the wall outlet. (Massage button is not depressed.)
- Use this appliance for its intended use as described in this manual.
- This unit is not a toy. Children should not use it or play with it.
- Never use while sleeping.
- Do not place or drop into water or any other liquid.

## Instructions for Use

To reduce the risk of electric shock, this appliance has a polarized plug (one blade is wider than the other). This plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install the proper outlet. Do not change the plug in any way.

1. It is highly recommended that Foot Pro Ultra is used with socked feet. Do not use the unit with shoes.
2. Place Foot Pro Ultra on the floor in front of your chair or seat. Make sure the unit is in the "OFF" position. Buttons are located on the bottom of the unit.

**Caution:**  
**Do not**  
**immerse in**  
**any liquid.**

3. Plug the unit into the power source (110-120V).
4. To activate massage action, press massage button. Control intensity between low and high by depressing Lo/Hi button located on left side of panel.
5. To activate infrared heat function, press red heat button on right side of unit. Allow a few moments for the heating action to warm up.
6. You can use the massage and heat functions together or separately.

This appliance is intended for household use only.

**NOTE:** Because the heat can be used without massage, make sure to turn the heat off after using the massager. Do not use the infrared heat continuously for more than 20 minutes. Do not use the infrared if you have sensitive skin, poor circulation, discolored areas, swelling, burning, inflammation or where skin eruptions or sores are present. Heat should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.

## Maintenance

**Caution:** All servicing of this unit must be performed by authorized HoMedics service personnel only. This appliance has no user serviceable parts.

- To Clean** Unplug the unit and allow it to cool before cleaning. Clean only with a soft, damp cloth. Never allow water or any other liquid to come into contact with the control switch.
- After cleaning, wipe surfaces with a dry cloth.
  - Never use abrasives, detergents or cleaning agents.
  - Periodically straighten cord if it becomes twisted.

**To Store** Unplug the appliance from the outlet and allow to cool.

- Store in its box or a clean, dry place.
- To avoid breakage, do not wrap the power cord around the unit.
- Wrap cord around the storage posts located on the bottom of the unit.
- Do not hang the unit by the power cord.

## Indulge your feet. Revitalize your body.

Thank you for purchasing Foot Pro Ultra, HoMedics' Luxury Foot Massager. Like all HoMedics products, it is built with high-quality craftsmanship and comes backed by our two-year limited warranty. Foot Pro Ultra is guaranteed to provide you with soothing massage relief for years to come.

HoMedics is dedicated to bringing you better health through innovation. However, you should never substitute this, or any other manufacturer's product, for professional medical attention. Persistent ailments and symptoms could indicate a more serious condition.

## Foot Pro Ultra Features

- Stimulating massage action rejuvenates the entire body
- Penetrating infrared heat helps soothe & relax tired feet
- Choose from gentle to invigorating massage action
- Convenient toe-tap massage control

# FootPro™ Ultra

## Luxury Foot Massager

### Massage Platform

With raised accunodules for all over stimulation

### Infrared Heat

